

Protein extraction, digestion, and TMT labeling

Lysis buffer:

RIPA lysis buffer kit, SC-24948. 10 mL
NP-40, 100 μ L
Protease inhibitor cocktail, 1 tablet
100 mM PMSF, 100 μ L

Cell lysis procedure

1. $\sim 10^7$ cells are transferred from cell culture dishes or flask to 1.5 or 2.0 mL Eppendorf tubes after washed with cold PBS x 2. Then, 200 μ L of lysis buffer solution is added to the cell pellets (Note: The volume of lysis buffer used varies depending on the amount of cells).
2. The cells are homogenized by sonication on ice by 30 strokes x 2.
3. The homogenized cells are incubated on ice by 1~2 hours.
4. After centrifugation at 14,000g for 10 min, the supernatant is pipetted to a clean 1.5 mL Eppendorf tube.
5. Another 200 μ L of lysis buffer solution is added to the cell pellets and steps 3 is repeated.
6. After centrifugation at 14,000g for 10 min, the supernatant is combined to the first supernatant.
7. Protein concentration is measured by BCA.

Protein digestion

1. 300 μ g of protein sample is diluted in 1mg/mL concentration by 50 mM Tris-HCl (pH 8.6) and 0.1% SDS buffer.
2. Proteins are reduced by incubation in 10 mM DTT at 50°C for 30 min followed by carboxymethylation with 25 mM iodoacetamide in the dark for 2 h.
3. Proteins are precipitated using cold acetone (85% or six volumes of sample volume) at - 20°C overnight. The protein pellet is generated by centrifugation at 15000 rpm for 10 min and the supernatant is removed using glass Pasteur pipette.
4. The protein pellet is dissolved in 500 μ l of 80 mM triethylammonium bicarbonate buffer. The protein solution should be clear after vortex or sonication. If not, add 1 or 2 μ g trypsin to digest for 4 hours. The protein concentration in each sample is measured again by BCA.
5. Divide the protein solution in 60~100 aliquots.
6. One aliquot of proteins (60 μ g) is digested using trypsin at a protein/enzyme ratio of 30:1 (by mass) at 37°C for overnights. Trypsin is added by two time-intervals, one (1 μ g) at the start and the other (1 μ g) after digestion for 4 hours.

TMT labeling

1. Immediately before use, equilibrate the TMT Label Reagents to room temperature. For the 0.8 mg vials, add 40 μ l of anhydrous acetonitrile to each tube. Allow the reagent to dissolve for 5 minutes with occasional vortexing. Briefly centrifuge the tube to gather the solution. **Note:** Reagents dissolved in anhydrous acetonitrile or DMSO are stable for one week when stored at -20°C and warmed to room temperature before opening.
2. Carefully add ~ 40 μ l of the TMT Label Reagent to each 25-100 μ g sample. Each sample is labeled with one specific TMT reagent (TMT-126 - TMT-131). Normally, duplicated samples are prepared and labeled with two different TMT reagents to have two analytical repeats. Alternatively, transfer the previously reduced and alkylated protein digests to the TMT Reagent vial. **Note:** A 100 μ l glass syringe or positive displacement pipette may be necessary to accurately measure and dispense TMT Reagents in volatile acetonitrile solvent.
3. Incubate the reaction for 2 hour at room temperature.
4. Add 15 μ l of 5% hydroxylamine to the sample and incubate for 15 minutes to quench the reaction.
5. Combine all of the TMT-labeled samples.
6. Desalt the combined sample using a C18/hypercarbon desalting column.
7. Run LC-MS/MS for the combined samples.

Alternatively with more cautions:

5. Desalt each sample using a C18/hypercarbon desalting column.
6. Run LC-MS/MS to determine the signal of β -actin for each sample.
7. Combine samples at equal amounts based on β -actin signals.
8. Run LC-MS/MS for the combined samples.

Note: For label-free proteomics, the TMT labeling step is skipped and the triethylammonium bicarbonate buffer is replaced by ammonium bicarbonate.